



Serological diagnosis of Lyme disease

Borrelia

- Gram-negative bacteria of the family Spirochaetaceae
- At least 36 genospecies are known, the group of pathogenic strains is referred to as *Borrelia burgdorferi sensu lato*

- Pathogenic species for humans:

Europe: *Borrelia burgdorferi sensu stricto*
Borrelia garinii
Borrelia afzelii
(*Borrelia spielmanii*)


USA: *Borrelia burgdorferi sensu stricto*

- Causative agent of Lyme disease (Borreliosis)




Ticks (genus Ixodes): Vectors of Lyme disease

- 3 developmental stages: larva, nymph, adult
- Blood feeding is necessary before every moult
- Transmission of *Borrelia* by ticks of all development stages, most frequently by nymphs
- Whole life cycle: \approx 3 years




Global distribution of ticks transmitting Borrelia



How do ticks become infected with Borrelia?

- Blood feeding on an already infected host
- Co-feeding in close vicinity to an infected tick
- Transovarial transmission to unfertilised eggs (\rightarrow infected larvae)
- Infection is kept for a lifetime



Transmission risk after tick bite (Germany)

- 743 ticks collected from humans (Maiwald et al., 1998)

57%	Nymphs
41%	Adult female Ticks
2%	Larvae

- 11% of the ticks were infected with *Borrelia*
- 27% of bites from infected ticks resulted in a human infection
- In total, \approx 4% of all tick bites lead to an infection



Infection prophylaxis (no vaccine available!):

1. Avoid tick exposure
2. After tick bite: avoid squeezing, remove tick quickly

Clinical symptoms and stages

Early local
Days to weeks

Disseminated
Weeks to months

Chronic
Months to years

Erythema migrans (EM)

Lymphocytoma

Facial paresis

Lyme arthritis

Acrodermatitis chronica atrophicans (ACA)

Problem: sometimes no EM, unspecific symptoms

Distribution of pathogenic Borrelia species

EUROPE	Ticks (n=90)	CSF (n=43)	Skin (n=68)	Synovial fluid (n=32)
B. burgdorferi s.s.	20%	19%	6%	33%
B. afzelii	9%	12%	84%	29%
B. garinii	71%	69%	10%	38%

Wilske et al. Diagnosis of Lyme borreliosis in Europe. Vector-Borne and Zoonotic Diseases.3(2003):215-227.

Laboratory diagnosis of Lyme disease

Direct detection methods: Microscopy, Culture, PCR
 ▶ low sensitivity (10-70%), time consuming, not standardised

Serological techniques: ELISA, Immunoblot, IFA

Stage	Symptoms	Sensitivity*	Comments
Local	Erythema migrans	20-50%	Predominantly IgM → Sensitivity can be significantly increased by detection of anti-VisE antibodies (IgG)
Disseminated	Neuroborreliosis	70-90%	IgM and IgG → with longer lasting disease predominance of IgG
Chronic	Lyme arthritis	90-100%	Normally only IgG

*taken from MIQ 12-2000

Serological two-step strategy (CDC, RKI, ...)

Screening IgG and IgM: ELISA or IFA (sensitive)

- positive or borderline → Confirmation
- negative → negative serological result

Confirmation IgG and IgM: Immunoblot (specific)

- positive → positive serological result
- borderline → doubtful serological result
- negative → negative serological result

VisE: the main antigen for Borrelia serology

VisE Variable major protein (VMP)-like sequence, Expressed

- 35 kDa surface lipoprotein of Borrelia only expressed *in vivo*!
- Constant VisE recombination protects the bacteria from elimination by the immune system
- Lyme disease patients show an early and strong IgG antibody response to VisE

Prevalences of IgG antibodies against VisE

Cohort	Anti-VisE ELISA IgG positive (%)
Culture-positive EM (n = 41)	63
Early disseminated (n = 12)	92
Acute neuroborreliosis (n = 17)	100
Lyme arthritis (n = 23)	87
Controls (n = 110)	2


Lawrenz et al., 1999 (group of Prof. Norris, University of Houston, Texas, USA)

Diagnostic value of VisE: ELISA

Erythema migrans n = 50		EUROIMMUN Anti-Borrelia plus VisE ELISA (IgG)		44%
		positive	negative	
full-extract ELISA without VisE	positive	22	0	70%
	negative	13	15	

⇒ VisE increases sensitivity in early stage by 26% !

Diagnostic value of VisE: Immunoblot



	B. afzelii IgG	B. afzelii plus VisE IgG	B. afzelii IgM	B. afzelii plus VisE IgM	Effect of VisE on blot sensitivity IgG + IgM
EM (n=47)	40%	62%	68%	70%	80 → 89%
NB (n=27)	78%	93%	48%	48%	85 → 96%
Arthritis (n=33)	94%	94%	15%	15%	94 → 94%
ACA (n=8)	100%	100%	13%	13%	100 → 100%

The addition of VisE leads to a significant increase of sensitivity!

Specificity of VisE

Disease	n	VisE IgG positive
Rheumatoid arthritis (RF seroneg.)	17	0 (0%)
Rheumatoid arthritis (RF seropos.)	23	1 (4%)
Finger polyarthrosis	25	0 (0%)
Fibrosyalgia	14	0 (0%)
Osteoarthritis	26	4 (15%)
Spondyloarthropathic arthritis (seroneg.)	22	1 (5%)
Total	119	6 (5%)


Specificity of anti-VisE antibodies investigated in patients with rheumatic diseases (Käfer et al., 2006)

Test system: Blot

Study group	Sera tested	Anti-VisE reactive sera (IgG)
Syphilis (secondary or latent)	24	0 (0%)
Louse-borne relapsing fever	11	1 (9%)
Oral infections	6	0 (0%)
Rheumatoid arthritis	7	0 (0%)
Healthy subjects	28	0 (0%)

Magnarelli et al., 2002

Test system: ELISA




European Patent Application: EP 2 199 303 A1

OspC advanced
A new designer antigen

OspC – Target of early IgM antibodies


OspC Outer surface protein C (p25)



- Most important target of IgM mediated immune response
- Upregulated during blood feeding
- Enables Bacteria to migrate into salivary glands


Starting point

native OspC (dimeric)




- Optimal antigen with high sensitivity and specificity
- Bivalent IgM binding with high affinity
→ low antigen input needed
- Very difficult to produce

recombinant OspC (monomeric)



- Easier production
- Monovalent IgM binding with low affinity
→ high antigen input needed
→ higher number of unspecific reactions

Objective: development of recombinant OspC with properties of native OspC




OspC advanced: Comparison with native OspC

Cohort	n	Prevalence of anti-OspC (IgM)					
		<i>B. burgdorferi</i>		<i>B. burgdorferi</i>		<i>B. garinii</i>	
		adv	native	adv	native	adv	native
Active borreliosis	188	85%	81%	88%	54%	60%	84%
Past infection (persisting IgM)	38	13%	13%	13%	13%	18%	19%
Acute EBV	30	0%	3%	0%	0%	3%	0%
Pregnant women	80	2%	2%	2%	2%	2%	2%
Blood donors	90	4%	4%	4%	4%	4%	2%

Scientific presentation at the 12th European Congress of Clinical Microbiology and Infectious Diseases (ECCMID), Milan, Italy, May 2010

Same sensitivity and specificity as native OspC!





OspC advanced: "Close to nature" rec. OspC

OspC advanced


- Combines advantages of native and recombinant OspC
- Same diagnostic characteristics as native OspC
- Easier to produce than native OspC → less costly, consistent quality
- 32% more specific than monomeric rec. OspC

N-terminal disulphide-bridging of Borrelia outer surface protein C increases its diagnostic and vaccine potentials

Christian Probst, Anthonina Ott, Thomas Scheper, Wolfgang Meyer, Winfried Stöcker, Lars Komorowski*



EUROPLUS (IgG, IgM): Anti-Borrelia IFA



Ig class specific EUROPLUS tests:

- Borrelia afzelii plus VisE antigen (IgG)
- Borrelia afzelii plus OspC antigen (IgM)



**Anti-Borrelia plus VisE (IgG)
Anti-Borrelia (IgM)**

Clinical studies with anti-Borrelia ELISAs

Cohort	n	Anti-Borrelia plus VisE IgG	Anti-Borrelia IgM	Anti-Borrelia IgG / IgM
Erythema migrans	205	76%	68%	91%
Neuroborreliosis	80	90%	49%	96%
Facial paresis	16	100%	50%	100%
Arthritis	49	84%	43%	94%
ACA	14	93%	21%	93%
Healthy blood donors	500	5%	2%	

Summarized results of various studies

EUROIMMUN Anti-Borrelia ELISA (IgG, IgM) agree 99% with the results from EQUALIS, INSTAND, IQS and Labquality quality assurance schemes

		Designated result from EQUALIS, INSTAND, IQS, Labquality (IgG)		
		positive	borderline	negative
EUROIMMUN Anti-Borrelia plus VisE ELISA (IgG)	positive	75	1	1
	borderline	0	2	1
	negative	0	0	51

		Designated result from EQUALIS, INSTAND, IQS, Labquality (IgM)		
		positive	borderline	negative
EUROIMMUN Anti-Borrelia ELISA (IgM)	positive	34	1	1
	borderline	0	1	3
	negative	0	1	110

Anti-Borrelia ELISAs in external QAS

EQUALIS: 10/2004-09/2010 | INSTAND: 10/2006-11/2010 | Labquality 02/2008-01/2011

	EQUALIS		INSTAND		Labquality							
	IgG	IgM	IgG	IgM	IgG	IgM						
EUROIMMUN	73	97	71	83	347	95	415	91	101	98	158	97
Genentech	191	96	166	88	767	84	893	97	106	98	152	88
Diagnostik LABOHY (CLIA)	125	94	125	93	36	100	41	80	194	94	179	94
Mikrogen	24	73	30	30	339	88	339	92	-	-	10	100
DAND / Oxoid	872	70	871	88	-	-	12	100	-	-	80	87
WilsonSerotec	-	-	-	-	255	84	388	92	-	-	-	-
Vitrocheck	-	-	-	-	333	85	341	93	43	93	25	88
Abnova	-	-	-	-	8	83	-	-	-	-	-	-
Bioss	-	-	-	-	-	-	-	-	43	87	45	87
IRL	-	-	-	-	-	-	-	-	21	83	21	80

n = number of participants | % = percentage of correct results



NEW

**Anti-Borrelia Select (IgG)
Anti-Borrelia Select (IgM)**

Tested Borrelia antigens

IgG

Antigens:

VisE
p17/p18
DbpA
OspC
BmpA (p39)
p100/p83

Species:

Borrelia burgdorferi
Borrelia garinii
Borrelia afzelii
Borrelia spielmanii

IgM

Antigens:

OspC – native purified
OspC – recombinant monomer
OspC – recombinant dimer (*advanced*)
VisE
BmpA (p39)
DbpA

Species:

Borrelia burgdorferi
Borrelia garinii
Borrelia afzelii
Borrelia spielmanii

Anti-Borrelia Select ELISA (IgG, IgM):

Mixture of recombinant antigens from various pathogenic *Borrelia* species

- IgG test with VisE
- IgM test with dimeric *OspC advanced* (optimised for use in ELISA)

Anti-Borrelia Select ELISA: Specificity

	Positive and grey zone results of ELISA			
	IgG		IgM	
	Lysate Ba, Bb, Bg plus VisE	Anti-Borrelia Select	Lysate Ba, Bb, Bg	Anti-Borrelia Select
Anti-Treponema positive (n=92)	65 (71 %)	3 (3 %)	26 (28 %)	6 (7 %)
Autoantibody positive (n=28)	6 (22 %)	5 (19 %)	7 (26 %)	1 (4 %)
Sera from daily laboratory routine (IgG: n=187 / IgM: n=185)	33 (18 %)	11 (6 %)	10 (5 %)	6 (3 %)

According to the literature, the prevalence of anti-Borrelia antibodies in healthy blood donors is 5 – 10% (IgG) and about 2% (IgM), respectively.

How about sensitivity?

Sera from external quality assessment INSTAND / 2004-2011

n=27	Anti-Borrelia Select (IgG)	Anti-Borrelia Select (IgM)
Sensitivity	100 %	100 %
Specificity	100 %	100 %

Sera from patients with erythema migrans (n=33)

Number of positive results		
Lysate Ba, Bb, Bg plus VisE (IgG)	27 (82 %)	30 (91 %)
Lysate Ba, Bb, Bg (IgM)	25 (76 %)	
Anti-Borrelia Select (IgG)	21 (64 %)	30 (91 %)
Anti-Borrelia Select (IgM)	22 (67 %)	

EUROIMMUN Anti-Borrelia ELISAs

Anti-Borrelia Screening ELISA (IgG + IgM)



Anti-Borrelia plus VisE (IgG)
Anti-Borrelia (IgM)

NEW


Anti-Borrelia Select (IgG)
Anti-Borrelia Select (IgM)

- Lysate based (*B. afzelii*, *B. burgdorferi*, *B. garinii*)
- IgG test with rec. **VisE**
- IgM test with high concentration of **native OspC**
- Full antigen spectrum
- Highest sensitivity

- Based on selection of highly specific recombinant *Borrelia* antigens
- IgG test with **VisE**
- IgM test with rec. **OspC advanced** (optimized for use in ELISA)
- Fewer cross reactivities (other infectious diseases, autoimmune diseases)


Confirmatory tests

EUROIMMUN immunoblots for confirmation



Anti-Borrelia EUROLINE-WB (IgG, IgM)

Lysate of *B. afzelii* plus membrane chip with VisE



- Combines advantages of native and recombinant antigens
- VisE as a sensitive, cross-species marker (IgG): already in the early phase of an infection
- Contains native OspC as main target of IgM antibodies
- Rapid at-a-glance evaluation:
Anti-Borrelia IgM: p25/OspC, Anti-Borrelia IgG: VisE
- Computer based evaluation and data management using the EUROlineScan software

Specificity of Anti-Borrelia EUROLINE-WB

Disease	n	Anti-Borrelia positive			VisE IgG positive
		IgG	IgM	IgG/IgM	
Rheumatoid arthritis (RF seropos.)	17	0 (0%)	0 (0%)	0 (0%)	0 (0%)
Rheumatoid arthritis (RF seroneg.)	23	1 (4%)	0 (0%)	1 (4%)	1 (4%)
Finger polyarthrosis	18	1 (6%)	0 (0%)	1 (6%)	0 (0%)
Fibromyalgia	14	0 (0%)	0 (0%)	0 (0%)	0 (0%)
Collagenosis	25	4 (16%)	0 (0%)	4 (16%)	4 (16%)
Spondyloarthropathy/psoriatic arthritis (seroneg.)	22	1 (5%)	0 (0%)	1 (5%)	1 (5%)
Total	118	7 (6%)	0 (0%)	7 (6%)	6 (5%)
Healthy blood donors	46	0 (12%)	0 (0%)	0 (12%)	0 (11%)


The EUROIMMUN Anti-Borrelia EUROLINE-WB (IgG, IgM) agree 100% with the results from the INSTAND quality assurance scheme

Sera investigated: 16 serologically precharacterized sera (IgG)
16 serologically precharacterized sera (IgM)
Origin of sera: INSTAND e. V., Düsseldorf, Germany


n = 16		Designated result from INSTAND (IgG)		
		positive	borderline	negative
EUROIMMUN Anti-Borrelia EUROLINE-WB (IgG)	positive	7	0	0
	borderline	0	0	0
	negative	0	0	0

n = 16		Designated result from INSTAND (IgM)		
		positive	borderline	negative
EUROIMMUN Anti-Borrelia EUROLINE-WB (IgM)	positive	2	0	0
	borderline	1	0	0
	negative	0	0	11

Lipid antigens



Anti-Borrelia EUROLINE-RN-AT



Performance of EUROLINE-RN-AT


Antigen	IgG	
	Prevalence*	Specificity*
VisE Ba	66%	99%
VisE Bb	89%	99%
VisE Bg	68%	95%
Lipid Ba	25%	100%
Lipid Bb	25%	100%
p83	54%	95%
p39	61%	99%
OspC	49%	96%
p58 (BB_A34)	21%	98%
p21 (BB_K53)	9%	99%
p20 (BB_Q03)	7%	100%
p19 (BB_N38)	9%	99%
p18 (BB_P38)	22%	99%

Evaluation of 617 sera from

- suspected LD
- clinically characterised LD
- control groups (blood donors, pregnant women, other infections)

Antigen	IgM	
	Prevalence*	Specificity*
VisE Bb	5%	99%
p39 Ba	16%	99%
OspC Ba	88%	99%
OspC Bb	77%	99%
OspC Bg	84%	97%

*referring to anti-Borrelia EUROLINE-WB



Anti-Borrelia CSF ELISA

Anti-Borrelia plus VisE ELISA (IgG) EI 2132-9601-L-G
Anti-Borrelia ELISA (IgM) EI 2132-9601-L-M

Calibration: 4 point / optional 5 (IgM) or 6 point (IgG)	Result interpretation:
Sample dilutions: Serum 1:404 / CSF 1:2	CSQrel. < 0.6 implausible result
Incubation: 60/60/15 (room temperature)	CSQrel. < 1.3 normal range
	CSQrel. 1.3 - 1.5 borderline range
	CSQrel. > 1.5 indication of pathogen-specific antibody production in CNS

- Excellent correlation of the calculated CSQrel. with clinically proved neuroborreliosis (sensitivity > 95 %)
- Over 95 % of patients with other neurological diseases have non-pathological Borrelia specific CSQrel. (specificity > 95 %)
- Best pass rates in the neuroborreliosis QAS by INSTAND*

**other manufacturers: Virion/Serion, Virotech, Mikrogen, Dade Behring, DiaSorin*

CSF diagnostics with Anti-Borrelia EL-RN-AT (IgG)

Interpretation of results

Depending on the number of specific bands (n) in CSF and serum:

- $n_{CSF} = n_{Serum}$ Indication of pathogen-specific antibody production in the CNS
- $n_{CSF} > n_{Serum}$ Interpretation based on the intensity of bands: (1) in CSF and serum possible
- $n_{CSF} = 0$ No indication of pathogen-specific antibody production in the CNS
- $n_{Serum} = 0$ Indication of pathogen-specific antibody production in the CNS

Advantages:


- High specificity (100%) and sensitivity (97%)
- Broad range of antigens
- Standardised dilution for CSF (1:4)
- Low sample volumes for CSF (250 µl)
- Short incubation times (≈ 300 min)
- Automated incubation (EBM) and evaluation (ELS)

Summary

- Lyme disease is a **clinical diagnosis**, based on symptoms!
- Serological **two-step strategy** is recommended: Sensitive screening + specific confirmation
- **VisE** (IgG antibodies) and **OspC** (IgM antibodies) are the most important antigens
- ELISA and Blot can be used with CSF samples to detect **intrathecally produced** antibodies
- **New EUROIMMUN tests:**
Anti-Borrelia Select ELISA (IgG, IgM)
Anti-Borrelia EUROLINE-RN-AT-adv (IgM)

Under Development

Antibody extraction from dried blood spots



Comparison blood spots ↔ serum

Cohort:


572 forest workers of various districts in Poland

	Dried blood spots vs. serum*	
	IgG ELISA	IgM ELISA
Sensitivity (%)	99.5	100
Specificity (%)	100	99.7
Positive predictive value (%)	100	99.1
Negative predictive value (%)	99.7	100


* excluding borderline results

Herbst et al., 12th International Conference on Lyme Borreliosis, Ljubljana, Slovenia (2010)

Antigen Source: Native Proteins



Antigen Source: Recombinant Proteins



Thank you!